

## Olerup SSP® HLA-A\*30

**Product number:** 101.429-12 – including *Taq* polymerase  
101.429-12u – without *Taq* polymerase  
**Lot number:** 2K3  
**Expiry date:** 2023-09-01  
**Number of tests:** 12  
**Number of wells per test:** 31+1

### CHANGES COMPARED TO THE PREVIOUS HLA-A\*30 LOT (1H7).

Well	5'-primer	3'-primer	rationale
7	Exchanged	-	5'-primer exchanged for the A*30:08:02 allele.

**THE NUMBER OF WELLS** is unchanged.

#### ALLELE COVERAGE:

All the HLA-A\*30 alleles, i.e. **A\*30:01 to A\*30:153 alleles**, recognized by the HLA Nomenclature Committee in April 2019<sup>1,2</sup> will be amplified by the primers in the HLA-A\*30 subtyping kit.

The HLA-A\*30 kit enables separation of the confirmed HLA-A\*30 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A\*30 kit also enables identification of null and alternatively expressed alleles.

The following HLA-A\*30 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
A*30:42, 30:81	26
A*30:56, 30:132N	27

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

#### RESOLUTION IN HLA-A\*30 HOMO- AND HETEROZYGOTES:

The A\*30:01,30:01, A\*30:01,30:02 and A\*30:02,30:02 genotypes give rise to unique amplification patterns.

#### INFLUENCE ON THE INTERPRETATION OF HLA-A\*30 SUBTYPINGS BY NON-HLA-A\*30 ALLELES:

None frequently occurring.

**HLA-A\*30**

**Release Note**

**Page 2 of 2**

**101.429-12 – including *Taq* polymerase**

**101.429-12u – without *Taq* polymerase**

**Lot No.: 2K3**

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.