

HLA-A\*31  
101.430-12 – licensed for PCR  
101.430-12u – not licensed for PCR  
Lot No.: Y30

1

## Release Note - Olerup SSP™ HLA-A\*31 – Lot No. Y30

Product number: 101.430-12 – licensed for PCR  
101.430-12u – not licensed for PCR  
Lot number: Y30  
Expiry date: 2009-September-01  
Number of tests: 12  
Number of tubes per test: 21

### CHANGES COMPARED TO THE PREVIOUS HLA-A\*31 LOT (V94):

Tube	5'-primer	3'-primer	rationale
14	New	New	New primer pair for A*3114N allele.
15	New	New	New primer pair for A*3115 allele.
16	New	New	New primer pair for A*3116 allele.
17	New	New	New primer pair for A*3117 allele.
18	New	New	New primer pair for A*3118 allele.
19	New	New	New primer pair for A*3119 allele.
20	New	New	New primer pair for A*3120 allele.
21	New	New	New primer pair for A*3121 allele.

THE NUMBER OF WELLS has been increased from 14 to 21.

#### ALLELE COVERAGE:

A\*3101 to A\*3121, i.e. all the currently recognized HLA-A\*31 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

#### RESOLUTION IN HLA-A\*31 HOMO- AND HETEROZYGOTES:

Very good.

#### INFLUENCE ON THE INTERPRETATION OF HLA-A\*31 SUBTYPINGS BY NON-HLA-A\*31 ALLELES:

None frequently occurring.

#### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.