

Olerup SSP[®] HLA-A*34

Product number: 101.425-06 – including *Taq* polymerase
101.425-06u – without *Taq* polymerase
Lot number: 23K
Expiry date: 2012-June-01
Number of tests: 6
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-A*34 LOT (58F):

Well	5'-primer	3'-primer	rationale
1	-	Modified	Primer modified for improved yield of specific PCR product.
2	-	Added	Primer added for improved yield of specific PCR product.
5	Modified	-	Improved specificity and yield of specific PCR product.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*34:01 to A*34:08, i.e. all the currently recognized HLA-A*34 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-A*34 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*34 SUBTYPINGS BY NON-HLA-A*34 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.