

Olerup SSP® HLA-A*34

Product number:	101.425-06 – including <i>Taq</i> polymerase 101.425-06u – without <i>Taq</i> polymerase
Lot number:	1H6
Expiry date:	2021-056-01
Number of tests:	6
Number of wells per test:	11+1

The HLA-A*34 primer set is unchanged compared to the previous *Olerup SSP*® HLA-A*34 (Lot No. 1F0).

THE NUMBER OF WELLS is unchanged

ALLELE COVERAGE:

A*34:01 to A*34:17, i.e. all the currently recognized HLA-A*34 alleles, will be amplified by the primers in the HLA-A*34 subtyping kit; www.ebi.ac.uk/imgt/hla, 2018-July-11¹, release 3.33.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-A*34 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*34 SUBTYPINGS BY NON-HLA-A*34 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.