

HLA-A*36
101.419-06 – licensed for PCR
101.419-06u – not licensed for PCR
Lot No.: Y23

1

Release Note - *Olerup* SSP™ HLA-A*36 – Lot No. Y23

Product number: 101.419-06 – licensed for PCR
101.419-06u – not licensed for PCR
Lot number: Y23
Expiry date: 2009-October-01
Number of tests: 6
Number of tubes per test: 6

CHANGES COMPARED TO THE PREVIOUS HLA-A*36 LOT (X04):

Tube	5'-primer	3'-primer	rationale
5	-	Modified	Increased resolution by avoiding co-amplification of the A*01 alleles.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*3601 to A*3604, i.e. all the currently recognized HLA-A*36 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN HLA-A*36 HOMO- AND HETEROZYGOTES:

All A*36 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-A*36 SUBTYPINGS BY NON-HLA-A*36 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.