

Olerup SSP[®] HLA-A*36

Product number:	101.419-06 – including <i>Taq</i> polymerase 101.419-06u – without <i>Taq</i> polymerase
Lot number:	89E
Expiry date:	2010-June-01
Number of tests:	6
Number of wells per test:	6

Changes compared to the previous HLA-A*36 Lot (Y23):

The HLA-A*36 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*3601 to A*3604, i.e. all the currently recognized HLA-A*36 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-A*36 HOMO- AND HETEROZYGOTES:

All A*36 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-A*36 SUBTYPINGS BY NON-HLA-A*36 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.