

## **Olerup SSP<sup>®</sup> HLA-A\*36**

**Product number:** 101.419-06 – including *Taq* polymerase  
101.419-06u – without *Taq* polymerase  
**Lot number:** 25G  
**Expiry date:** 2011-October-01  
**Number of tests:** 6  
**Number of wells per test:** 6

### **Changes compared to the previous HLA-A\*36 Lot (89E):**

The HLA-A\*36 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

A\*3601 to A\*3604, i.e. all the currently recognized HLA-A\*36 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-July-17, release 2.26.0.

### **RESOLUTION IN HLA-A\*36 HOMO- AND HETEROZYGOTES:**

All A\*36 genotypes give rise to unique amplification patterns.

### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*36 SUBTYPINGS BY NON-HLA-A\*36 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.