

Olerup SSP[®] HLA-A*36

Product number: 101.419-06 – including *Taq* polymerase
101.419-06u – without *Taq* polymerase
Lot number: 98K
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 7

CHANGES COMPARED TO THE PREVIOUS HLA-A*36 LOT (25G):

Well	5'-primer	3'-primer	rationale
7	New	New	New primer pair for the A*36:05 allele.

THE NUMBER OF WELLS has been increased from 6 to 7.

ALLELE COVERAGE:

A*36:01 to A*36:05, i.e. all the currently recognized HLA-A*36 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-A*36 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*36 SUBTYPINGS BY NON-HLA-A*36 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.