

HLA-A*66
101.427-06 – licensed for PCR
101.427-06u – not licensed for PCR
Lot No.: Y03

1

Release Note - Olerup SSP™ HLA-A*66 – Lot No. Y03

Product number: 101.427-06 – licensed for PCR
101.427-06u – not licensed for PCR
Lot number: Y03
Expiry date: 2009-July-01
Number of tests: 6
Number of tubes per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-A*66 LOT (V91):

Tube	5'-primer	3'-primer	rationale
2	Modified	-	Primer modified for increased specificity.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*6601 to A*6606, i.e. all the currently recognized HLA-A*66 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN HLA-A*66 HOMO- AND HETEROZYGOTES:

Good, except that that the A*6601,6602 genotype does not give rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-A*66 SUBTYPINGS BY NON-HLA-A*66 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.