

Lot No.: **08L**

Olerup SSP[®] HLA-A*66

Product number: 101.427-06 – including *Taq* polymerase
101.427-06u – without *Taq* polymerase
Lot number: 08L
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 14

CHANGES COMPARED TO THE PREVIOUS HLA-A*66 LOT (34G):

Well	5'-primer	3'-primer	rationale
3	-	Added	Primer added for the A*66:15 allele.
9	-	Added	Primer added for the A*66:14 allele.
11	-	Added	Primer added for the A*66:14 allele.
13	New	New	New primer pairs for the A*66:11 and A*66:13 alleles.
14	New	New	New primer pair for the A*66:12 allele.

THE NUMBER OF WELLS has been increased from 12 to 14 wells.

ALLELE COVERAGE:

A*66:01 to A*66:15, i.e. all the currently recognized HLA-A*66 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-A*66 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*66 SUBTYPINGS BY NON-HLA-A*66 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.