

101.844-12 – including *Taq* polymerase101.844-12u – without *Taq* polymerase

Lot No.: 30Y

Olerup SSP® HLA-A*68 Add-on

Product number: 101.844-12 – including *Taq* polymerase
101.844-12u – without *Taq* polymerase

Lot number: 30Y

Expiry date: 2017-October-01

Number of tests: 12

Number of wells per test: 2+1

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*68:01 and 68:11N alleles give different patterns in the HLA-A*68 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2015-January-19, release 3.19.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*68 ADD-ON SUBTYPINGS BY NON-HLA-A*68 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.