

101.844-12 – including *Taq* polymerase101.844-12u – without *Taq* polymeraseLot No.: **2E3****Olerup SSP® HLA-A*68 Add-on**

Product number:	101.844-12 – including <i>Taq</i> polymerase 101.844-12u – without <i>Taq</i> polymerase
Lot number:	2E3
Expiry date:	2019-02-01
Number of tests:	12
Number of wells per test:	2+1

The HLA-A*68 Add-on primer set is unchanged compared to the previous *Olerup SSP® HLA-A*68 Add-on (Lot No. 30Y)*.

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*68:01 and 68:11N alleles give different patterns in the HLA-A*68 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2016-April-15, release 3.24.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*68 ADD-ON SUBTYPINGS BY NON-HLA-A*68 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.