

Olerup SSP[®] HLA-A*80

Product number:	101.434-06 – including <i>Taq</i> polymerase 101.434-06u – without <i>Taq</i> polymerase
Lot number:	06K
Expiry date:	2012-May-01
Number of tests:	6
Number of wells per test:	3

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

A*80:01 and A*80:02, i.e. all the currently recognized HLA-A*80 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-A*80 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*80 SUBTYPINGS BY NON-HLA-A*80 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.