

Olerup SSP[®] HLA-A*80

Product number: 101.434-06 – including *Taq* polymerase
101.434-06u – without *Taq* polymerase
Lot number: 70M
Expiry date: 2014-March-01
Number of tests: 6
Number of wells per test: 3

CHANGES COMPARED TO THE PREVIOUS HLA-A*80 LOT (06K):

The HLA-A*80 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*80:01 and A*80:02, i.e. all the currently recognized HLA-A*80 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-July-14, release 3.5.0.

RESOLUTION IN HLA-A*80 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*80 SUBTYPINGS BY NON-HLA-A*80 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.