

101.513-24/04 – including *Taq* polymerase
101.513-24u/04u – without *Taq* polymerase

Lot No.: **24E**

Olerup SSP[®] HLA-B*08

Product number: 101.513-24/04 – including *Taq* pol.
101.513-24u/04u - without *Taq* pol.
Lot number: 24E
Expiry date: 2010-February-01
Number of tests: 24 tests – Product No. 101.513-24
4 tests – Product No. 101.513-04
Number of wells per test: 32

CHANGES COMPARED TO THE PREVIOUS HLA-B*08 LOT (X33):

Well	5'-primer	3'-primer	rationale
8	Modified	-	Increased specificity of specific primer pair.
29	New	New	New primer pair for the B*0832 allele.
30	New	New	New primer pair for the B*0833 allele.
31	New	New	New primer pair for the B*0834 allele.
32	New	New	New primer pair for the B*0835 allele.

THE NUMBER OF WELLS has been increased from 28 to 32.

ALLELE COVERAGE:

B*0801 to B*0835, i.e. all the currently recognized HLA-B*08 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN HLA-B*08 HOMO- AND HETEROZYGOTES:

The B*0801,0801 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*08 SUBTYPINGS BY NON-HLA-B*08 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.