

Olerup SSP[®] HLA-B*13

Product number: 101.515-12 – including *Taq* polymerase
101.515-12u - without *Taq* polymerase
Lot number: 97F
Expiry date: 2011-July-01
Number of tests: 12
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-B*13 LOT (96E):

Well	5'-primer	3'-primer	rationale
5	-	Added	Primer added for the B*1322 and B*1324 alleles.
6	Modified	-	Modified 5'-primer for increased yield of specific PCR product.
15	Added	Added	Primer pair added for the B*1322 and B*1323 alleles, and exchanged positive control primer pair.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*1301 to B*1324, i.e. all the currently recognized HLA-B*13 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-May-11, release 2.25.2.

RESOLUTION IN HLA-B*13 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*13 SUBTYPINGS BY NON-HLA-B*13 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.