

Olerup SSP® HLA-B*15 Add-on

Product number:	101.854-12 – including <i>Taq</i> polymerase 101.854-12u – without <i>Taq</i> polymerase
Lot number:	8E5
Expiry date:	2019-07-01
Number of tests:	12
Number of wells per test:	2+1

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-B*15:01:01:02N and 15:01:01G alleles give different patterns in the HLA-B*15 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2016-July-14, release 3.25.0.

The HLA-B*15 Add-on kit cannot distinguish these silent mutations: the B*15:01:01:01, 15:01:01:03-15:01:04 and 15:01:06-15:01:39 alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF B*15 ADD-ON SUBTYPINGS BY NON-HLA-B*15 ALLELES:

None of importance.