

Release Note - *Olerup* SSP™ HLA-B*27 – Lot No. **X94**

Product number: 101.521-24/04 – licensed for PCR
101.521-24u/04u - not licensed for PCR
Lot number: X94
Expiry date: 2009-June-01
Number of tests: 24 test – Product No. 101.521-24
4 tests – Product No. 101.521-04
Number of tubes per test: 24

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 LOT (V40):

| Tube | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|---|
| 20 | - | Added | Primer added for the B*2736 allele. |
| 22 | - | Added | Primer added for the B*2737 allele. |
| 24 | Added | Added | Primer pair added for separating the B*2729 and B*3702 alleles. |

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*2701 to B*2738, i.e. all the currently recognized HLA-B*27 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN HLA-B*27 HOMO- AND HETEROZYGOTES:

The B*2705,2705 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 SUBTYPINGS BY NON-HLA-B*27 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.