

101.565-12 – including *Taq* polymerase101.565-12u – without *Taq* polymeraseLot No.: **80F****Olerup SSP<sup>®</sup> HLA-B\*38**

**Product number:** 101.565-12 – including *Taq* polymerase  
101.565-12u – without *Taq* polymerase  
**Lot number:** 80F  
**Expiry date:** 2011-April-01  
**Number of tests:** 12  
**Number of wells per test:** 16

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*38 LOT (45E):**

Well	5'-primer	3'-primer	rationale
12	Modified	-	Increased yield of specific PCR product.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*3801 to B\*3817 i.e. all the currently recognized HLA-B\*38 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-January-16, release 2.24.0.

**RESOLUTION IN HLA-B\*38 HOMO- AND HETEROZYGOTES:**

The B\*3801,3801, B\*3801,3802 and B\*3802,3802 genotypes give rise to unique amplification patterns.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*38 SUBTYPINGS BY NON-HLA-B\*38 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.