

Olerup SSP[®] HLA-B*39 Add-on

Product number:	101.853-12 – including <i>Taq</i> polymerase 101.853-12u – without <i>Taq</i> polymerase
Lot number:	24N
Expiry date:	2014-August-01
Number of tests:	12
Number of wells per test:	2

THE NUMBER OF WELLS is 2 in an 8 well plate.

ALLELE COVERAGE:

The HLA-B*39:01 and 39:01:02:01L alleles give different patterns in the HLA-B*39 Add-on kit²; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

¹The HLA-B*39 Add-on kit cannot distinguish the silent mutations in the B*39:01:01:01, 39:01:03-39:01:08 and 39:01:10-39:01:12 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF B*39 ADD-ON SUBTYPINGS BY NON-HLA-B*39 ALLELES:

None of importance.