

## **Olerup SSP® HLA-B\*39 Add-on**

Product number:	101.853-12 – including <i>Taq</i> polymerase 101.853-12u – without <i>Taq</i> polymerase
Lot number:	46S
Expiry date:	2015-December-01
Number of tests:	12
Number of wells per test:	2

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*39 ADD-ON LOT (24N):**

The HLA-C\*07 Add-on primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is 2 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-B\*39:01 and 39:01:02:01L alleles give different patterns in the HLA-B\*39 Add-on kit<sup>2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-April-17, release 3.12.0.

<sup>1</sup>The HLA-B\*39 Add-on kit cannot distinguish the silent mutations in the B\*39:01:01:01, 39:01:01:03, 39:01:03-39:01:08, 39:01:10-39:01:14 and 39:01:16 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF B\*39 ADD-ON SUBTYPINGS BY NON-HLA-B\*39 ALLELES:**

None of importance.