

**101.853-12 – including *Taq* polymerase**

**101.853-12u – without *Taq* polymerase**

**Lot No.: 8G0**

## **Olerup SSP® HLA-B\*39 Add-on**

<b>Product number:</b>	<b>101.853-12 – including <i>Taq</i> polymerase</b> <b>101.853-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>8G0</b>
<b>Expiry date:</b>	<b>2022-09-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>2+1</b>

The HLA-B\*39 Add-on primer set is unchanged compared to the previous *Olerup SSP® HLA-B\*39 Add-on (Lot No. 4E2)*.

**THE NUMBER OF WELLS** is 3 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-B\*39:01 and 39:01:02:01L alleles give different patterns in the HLA-B\*39 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2018-April-16, release 3.32.0.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF B\*39 ADD-ON SUBTYPINGS BY NON-HLA-B\*39**

#### **ALLELES:**

None of importance.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.