

***Olerup SSP*[®] HLA-B*42**

Product number:	101.543-06 – including <i>Taq</i> polymerase 101.543-06u – without <i>Taq</i> polymerase
Lot number:	51E
Expiry date:	2010-May-01
Number of tests:	6
Number of wells per test:	10

CHANGES COMPARED TO THE PREVIOUS HLA-B*42 LOT (X78):

The HLA-B*42 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*4201 to B*4209, i.e. all the currently recognized HLA-B*42 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-B*42 HOMO- AND HETEROZYGOTES:

The B*4201,4201 and B*4201,4208 genotypes give rise to identical amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-B*42 SUBTYPINGS BY NON-HLA-B*42 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.