

HLA-B*45
101.564-06 – including *Taq* polymerase
101.564-06u – without *Taq* polymerase
Lot No.: **16E**

Release Note

Page 1 of 1

Olerup SSP[®] HLA-B*45

Product number:	101.564-06 – including <i>Taq</i> polymerase 101.564-06u – without <i>Taq</i> polymerase
Lot number:	16E
Expiry date:	2010-February-01
Number of tests:	6
Number of wells per test:	8

CHANGES COMPARED TO THE PREVIOUS HLA-B*45 LOT (X50):

The HLA-B*45 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*4501 to B*4507, i.e. all the currently recognized HLA-B*45 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN HLA-B*45 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*45 SUBTYPINGS BY NON-HLA-B*45 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.