

Olerup SSP[®] HLA-B*45

Product number: 101.564-06 – including *Taq* polymerase
101.564-06u – without *Taq* polymerase
Lot number: 86F
Expiry date: 2011-May-01
Number of tests: 6
Number of wells per test: 10

CHANGES COMPARED TO THE PREVIOUS HLA-B*45 LOT (16E):

Well	5'-primer	3'-primer	rationale
6	-	-	Exchanged positive control primer pair.
9	New	New	New primer pair for the B*4508 allele.
10	New	New	New primer pair for the B*4509 allele.

THE NUMBER OF WELLS has been increased from 8 to 10.

ALLELE COVERAGE:

B*4501 to B*4509, i.e. all the currently recognized HLA-B*45 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-May-11, release 2.25.2.

RESOLUTION IN HLA-B*45 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*45 SUBTYPINGS BY NON-HLA-B*45 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.