

101.564-06 – including *Taq* polymerase101.564-06u – without *Taq* polymeraseLot No.: **61M****Olerup SSP<sup>®</sup> HLA-B\*45**

**Product number:** 101.564-06 – including *Taq* polymerase  
101.564-06u – without *Taq* polymerase  
**Lot number:** 61M  
**Expiry date:** 2014-March-01  
**Number of tests:** 6  
**Number of wells per test:** 14

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*45 LOT (24K):**

Well	5'-primer	3'-primer	rationale
13	New	New	New primer pair for the B*45:12 allele.
14	New	New	New primer pair for the B*45:13 allele.

**THE NUMBER OF WELLS** has been increased from 12 to 14.

**ALLELE COVERAGE:**

B\*45:01 to B\*45:13, i.e. all the currently recognized HLA-B\*45 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-July-14, release 3.5.0.

**RESOLUTION IN HLA-B\*45 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*45 SUBTYPINGS BY NON-HLA-B\*45 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.