

Olerup SSP® HLA-B*45

Product number:	101.564-06 – including <i>Taq</i> polymerase 101.564-06u – without <i>Taq</i> polymerase
Lot number:	0G4
Expiry date:	2020-08-01
Number of tests:	6
Number of wells per test:	15+1

The HLA-B*45 primer set is unchanged compared to the previous *Olerup SSP*® HLA-B*45 (Lot No. 7D5).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*45:01 to B*45:21, i.e. all the currently recognized HLA-B*45 alleles, will be amplified by the primers in the HLA-B*45 SSP kit; www.ebi.ac.uk/imgt/hla, 2017-October-27, release 3.30.0.

The HLA-B*45 kit enables separation of the confirmed HLA-B*45 alleles as listed in the IMGT/HLA database 3.23.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B*45 kit also enables identification of null and alternatively expressed alleles.

RESOLUTION IN HLA-B*45 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*45 SUBTYPINGS BY NON-HLA-B*45 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.