

101.544-06 – including *Taq* polymerase101. 544-06u – without *Taq* polymeraseLot No.: **93M****Olerup SSP[®] HLA-B*46**

Product number: 101.544-06 – including *Taq* polymerase
101.544-06u – without *Taq* polymerase
Lot number: 93M
Expiry date: 2014-May-01
Number of tests: 6
Number of wells per test: 24

CHANGES COMPARED TO THE PREVIOUS HLA-B*46 LOT (29M):

Well	5'-primer	3'-primer	rationale
23	New	New	New primer pair for the B*46:28 allele.
24	New	New	New primer pair for the B*46:30 allele.

THE NUMBER OF WELLS is increased from 22 to 24.

ALLELE COVERAGE:

B*46:01 to B*46:30, i.e. all the currently recognized HLA-B*46 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-October-10, release 3.6.0.

RESOLUTION IN HLA-B*46 HOMO- AND HETEROZYGOTES:

The B*46:01,46:01 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*46 SUBTYPINGS BY NON-HLA-B*46 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.