

HLA-B\*47  
101.545-06 – including *Taq* polymerase  
101.545-06u – without *Taq* polymerase  
Lot No.: **13E**

Release Note

Page 1 of 1

## **Olerup SSP<sup>®</sup> HLA-B\*47**

<b>Product number:</b>	<b>101.545-06 – including <i>Taq</i> polymerase</b> <b>101.545-06u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>13E</b>
<b>Expiry date:</b>	<b>2010-February-01</b>
<b>Number of tests:</b>	<b>6</b>
<b>Number of Wells per test:</b>	<b>7</b>

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*47 LOT (X06):**

The HLA-B\*47 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

B\*4701 to B\*4705, i.e. all the currently recognized HLA-B\*47 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-January-11, release 2.20.0.

### **RESOLUTION IN HLA-B\*47 HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*47 SUBTYPINGS BY NON-HLA-B\*47 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.