

Lot No.: **10K**

Olerup SSP[®] HLA-B*47

Product number: 101.545-06 – including *Taq* polymerase
101.545-06u – without *Taq* polymerase
Lot number: 10K
Expiry date: 2012-May-01
Number of tests: 6
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-B*47 LOT (84F):

Well	5'-primer	3'-primer	rationale
5	Added	Added	Primer pair added for the B*47:06 allele, exchanged positive control primer pair
8	New	New	New primer pair for the B*47:07 allele.

THE NUMBER OF WELLS has been increased from 7 to 8.

ALLELE COVERAGE:

B*47:01 to B*47:07, i.e. all the currently recognized HLA-B*47 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-B*47 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*47 SUBTYPINGS BY NON-HLA-B*47 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.