

HLA-B*49
101.547-06 – including *Taq* polymerase
101.547-06u – without *Taq* polymerase
Lot No.: **17E**

Release Note

Page 1 of 1

Olerup SSP[®] HLA-B*49

Product number:	101.547-06 – including <i>Taq</i> polymerase 101.547-06u – without <i>Taq</i> polymerase
Lot number:	17E
Expiry date:	2010-February-01
Number of tests:	6
Number of wells per test:	8

CHANGES COMPARED TO THE PREVIOUS HLA-B*49 LOT (X57):

The HLA-B*49 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*4901 to B*4905, i.e. all the currently recognized HLA-B*49 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN HLA-B*49 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*49 SUBTYPINGS BY NON-HLA-B*49 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.