

Lot No.: **64F**

Olerup SSP[®] HLA-B*49

Product number: 101.547-06 – including *Taq* polymerase
101.547-06u – without *Taq* polymerase
Lot number: 64F
Expiry date: 2011-February-01
Number of tests: 6
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-B*49 LOT (17E):

Well	5'-primer	3'-primer	rationale
2	Modified	-	Increases specificity of specific primer pair.
4	Modified	-	Increases specificity of specific primer pair.
8	New	New	New primer pair for the B*4905 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*4901 to B*4905, i.e. all the currently recognized HLA-B*49 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-January-16, release 2.24.0.

RESOLUTION IN HLA-B*49 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*49 SUBTYPINGS BY NON-HLA-B*49 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.