

101.547-06 – including *Taq* polymerase101.547-06u – without *Taq* polymeraseLot No.: **05H****Olerup SSP<sup>®</sup> HLA-B\*49**

**Product number:** 101.547-06 – including *Taq* polymerase  
101.547-06u – without *Taq* polymerase  
**Lot number:** 05H  
**Expiry date:** 2012-March-01  
**Number of tests:** 6  
**Number of wells per test:** 11

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*49 LOT (64F):**

Well	5'-primer	3'-primer	rationale
2	Added	Added	Primer pair added for the B*4906 allele.
4	Added	Added	Primer pair added for the B*4909 allele.
9	New	New	New primer pair for the B*4907 allele.
10	New	New	New primer pair for the B*4908 allele.
11	New	New	New primer pair for improved allelic resolution.

**THE NUMBER OF WELLS** is increased from 8 to 11.

**ALLELE COVERAGE:**

B\*4901 to B\*4909, i.e. all the currently recognized HLA-B\*49 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-January-15, release 2.28.0.

**RESOLUTION IN HLA-B\*49 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*49 SUBTYPINGS BY NON-HLA-B\*49 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.