

101.548-06 – including *Taq* polymerase101.548-06u – without *Taq* polymeraseLot No.: **89F****Olerup SSP® HLA-B\*50**

**Product number:** 101.548-06 – including *Taq* polymerase  
101.548-06u – without *Taq* polymerase  
**Lot number:** 89F  
**Expiry date:** 2011-May-01  
**Number of tests:** 6  
**Number of wells per test:** 7

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*50 LOT (33E):**

Well	5'-primer	3'-primer	rationale
7	New	New	New primer pair for improved allelic resolution.

**THE NUMBER OF WELLS** has been increased from 6 to 7.

**ALLELE COVERAGE:**

B\*5001 to B\*5004, i.e. all the currently recognized HLA-B\*50 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-May-11, release 2.25.2.

**RESOLUTION IN HLA-B\*50 HOMO- AND HETEROZYGOTES:**

Very good.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*50 SUBTYPINGS BY NON-HLA-B\*50 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.