

Olerup SSP[®] HLA-B*50

Product number: 101.548-06 – including *Taq* polymerase
101.548-06u – without *Taq* polymerase
Lot number: 67M
Expiry date: 2014-March-01
Number of tests: 6
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-B*50 LOT (49K):

Well	5'-primer	3'-primer	rationale
3	-	-	Exchanged positive control primer pair.
13	New	New	New primer pairs for the B*50:10 and 50:12 alleles.
14	New	New	New primer pair for the B*50:11 allele.
15	New	New	New primer pair for the B*50:13 allele.
16	New	New	New primer pair for the B*50:14 allele.

THE NUMBER OF WELLS has been increased from 12 to 16.

ALLELE COVERAGE:

B*50:01 to B*50:14, i.e. all the currently recognized HLA-B*50 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-July-14, release 3.5.0.

RESOLUTION IN HLA-B*50 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*50 SUBTYPINGS BY NON-HLA-B*50 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.