

## **Olerup SSP<sup>®</sup> HLA-B\*50**

**Product number:** 101.548-06 – including *Taq* polymerase  
101.548-06u – without *Taq* polymerase  
**Lot number:** 33E  
**Expiry date:** 2010-March-01  
**Number of tests:** 6  
**Number of wells per test:** 6

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*50 LOT (X58):**

The HLA-B\*50 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

B\*5001 to B\*5004, i.e. all the currently recognized HLA-B\*50 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-January-11, release 2.20.0.

### **RESOLUTION IN HLA-B\*50 HOMO- AND HETEROZYGOTES:**

Very good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*50 SUBTYPINGS BY NON-HLA-B\*50 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.