

101.851-12 – including *Taq* polymerase

101.851-12u – without *Taq* polymerase

Lot No.: **35V**

Olerup SSP® HLA-B*51:11N

Product number: 101.851-12 – including *Taq* polymerase
101.851-12u – without *Taq* polymerase

Lot number: 35V

Expiry date: 2016-August-01

Number of tests: 12

Number of wells per test: 2+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*51:11N LOT (96R):

Well	5'-primer	3'-primer	rationale
3	-	-	Negative Control

THE NUMBER OF WELLS is increased from 2 to 3 wells

ALLELE COVERAGE:

The HLA-B*51:11N allele will give rise to a unique amplification pattern by the primers in the HLA-B*51:11N kit; www.ebi.ac.uk/imgt/hla, 2013-October-11, release 3.14.0.

RESOLUTION IN HLA-B*51:11N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*51:11N SUBTYPINGS BY NON-HLA-B*51:11N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.