

101.851-12 – including *Taq* polymerase

101.851-12u – without *Taq* polymerase

Lot No.: **0D7**

## **Olerup SSP® HLA-B\*51:11N**

**Product number:** 101.851-12 – including *Taq* polymerase  
101.851-12u – without *Taq* polymerase

**Lot number:** 0D7

**Expiry date:** 2018-April-01

**Number of tests:** 12

**Number of wells per test:** 2+1

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*51:11N LOT (35V):**

Well	5'-primer	3'-primer	rationale
3	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The HLA-B\*51:11N allele will give rise to a unique amplification pattern by the primers in the HLA-B\*51:11N kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-August-11, release 3.21.1.

### **RESOLUTION IN HLA-B\*51:11N HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*51:11N SUBTYPINGS BY NON-HLA-B\*51:11N ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.