

HLA-B*51:11N

101.851-12 – including *Taq* polymerase101.851-12u – without *Taq* polymeraseLot No.: **4K5****Olerup SSP[®] HLA-B*51:11N**

Product number:	101.851-12 – including <i>Taq</i> polymerase 101.851-12u – without <i>Taq</i> polymerase
Lot number:	4K5
Expiry date:	2023-10-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*51:11N LOT (4F3):

The HLA-B*51:11N primer set is unchanged compared to the previous *Olerup SSP[®]* HLA-B*51:11N (Lot No. 4F3).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-B*51:11N allele will give rise to a unique amplification pattern by the primers in the HLA-B*51:11N kit; www.ebi.ac.uk/imgt/hla, 2019-July-10, release 3.37.0.

RESOLUTION IN HLA-B*51:11N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*51:11N SUBTYPINGS BY NON-HLA-B*51:11N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.