

Olerup SSP[®] HLA-B*52

Product number: 101.562-06 – including *Taq* polymerase
101.562-06u – without *Taq* polymerase
Lot number: 15G
Expiry date: 2011-September-01
Number of tests: 6
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-B*52 LOT (90E):

Well	5'-primer	3'-primer	rationale
4	Modified	-	Modified primer for increased specificity.
9	Added	Added	Added primer pair for the B*5214 allele.
16	New	New	New primer pair for the B*5213 allele.

THE NUMBER OF WELLS has been increased from 15 to 16.

ALLELE COVERAGE:

B*5201 to B*5214, i.e. all the currently recognized HLA-B*52 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HLA-B*52 HOMO- AND HETEROZYGOTES:

The B*5201,5201 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*52 SUBTYPINGS BY NON-HLA-B*52 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.