

## Release Note - Olerup SSP™ HLA-B\*52 – Lot No. Y31

Product number: 101.562-06 – licensed for PCR  
101.562-06u – not licensed for PCR  
Lot number: Y31  
Expiry date: 2009-August-01  
Number of tests: 6  
Number of tubes per test: 14

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*52 LOT (X01):

Tube	5'-primer	3'-primer	rationale
2	-	Exchanged	More specific B*5209 primer pair.
3	Modified	-	Increased specificity of specific primer pair
4	Exchanged	Exchanged	B*5202-specific primer pair.
5	Exchanged	-	More specific primer pair.
6	Exchanged	-	More specific primer pair-
9	Exchanged	Exchanged	B*520601-specific primer pair.
10	-	-	Previously primer pair 15.
15	Moved	Moved	Moved to vial 10.

THE NUMBER OF WELLS has been decreased from 15 to 14.

### ALLELE COVERAGE:

B\*5201 to B\*5211, i.e. all the currently recognized HLA-B\*52 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

### RESOLUTION IN HLA-B\*52 HOMO- AND HETEROZYGOTES:

The B\*5201,5201 genotype gives rise to a unique amplification pattern.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*52 SUBTYPINGS BY NON-HLA-B\*52 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.