

Olerup SSP[®] HLA-B*54

Product number: 101.569-06 – including *Taq* polymerase
101.569-06u – without *Taq* polymerase
Lot number: 41F
Expiry date: 2010-December-01
Number of tests: 6
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-B*54 LOT (Y26):

Well	5'-primer	3'-primer	rationale
15	New	New	New primer pair for the B*5414, B*5415 and B*5416 alleles.
16	New	New	New primer pair for the B*5417 allele.

THE NUMBER OF WELLS has been increased from 14 to 16.

ALLELE COVERAGE:

B*5401 to B*5417, i.e. all the currently recognized HLA-B*54 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008 October-10, release 2.23.0.

RESOLUTION IN HLA-B*54 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*54 SUBTYPINGS BY NON-HLA-B*54 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.