

101.569-06 – including *Taq* polymerase101.569-06u – without *Taq* polymeraseLot No.: **08H****Olerup SSP® HLA-B\*54**

**Product number:** 101.569-06 – including *Taq* polymerase  
101.569-06u – without *Taq* polymerase  
**Lot number:** 08H  
**Expiry date:** 2012-April-01  
**Number of tests:** 6  
**Number of wells per test:** 16

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*54 LOT (41F):**

Well	5'-primer	3'-primer	rationale
4	Added	-	Increased yield of specific PCR product, primer added for the B*5419 allele.
11	Added	Added	Primer pair added for the B*5418 allele.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*5401 to B\*5417, i.e. all the currently recognized HLA-B\*54 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-January-16, release 2.28.0.

**RESOLUTION IN HLA-B\*54 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*54 SUBTYPINGS BY NON-HLA-B\*54 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.