

101.569-06 – including *Taq* polymerase101.569-06u – without *Taq* polymeraseLot No.: **30M****Olerup SSP<sup>®</sup> HLA-B\*54**

**Product number:** 101.569-06 – including *Taq* polymerase  
101.569-06u – without *Taq* polymerase  
**Lot number:** 30M  
**Expiry date:** 2013-December-01  
**Number of tests:** 6  
**Number of wells per test:** 20

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*54 LOT (08H):**

Well	5'-primer	3'-primer	rationale
3	-	Modified	Improved specificity.
11	Moved	Moved	Primer pair moved to well 16, for resolution of the B*54:18 and B*54:22 alleles.
16	Added	Added	Primer pair from well 11.
17	New	New	New primer pair for the B*54:20 allele.
18	New	New	New primer pair for the B*54:21 allele.
19	New	New	New primer pair for the B*54:23 allele.
20	New	New	New primer pair for the B*54:24 allele.

**THE NUMBER OF WELLS** is increased from 16 to 20.

**ALLELE COVERAGE:**

B\*54:01 to B\*54:24, i.e. all the currently recognized HLA-B\*54 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-April-08, release 3.4.0.

**RESOLUTION IN HLA-B\*54 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*54 SUBTYPINGS BY NON-HLA-B\*54 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.