

Olerup SSP[®] HLA-B*57

Product number: 101.567-12 – including *Taq* polymerase
101.567-12u - without *Taq* polymerase
Lot number: 59F
Expiry date: 2011-January-01
Number of tests: 12
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-B*57 LOT (62E):

Well	5'-primer	3'-primer	rationale
13	Exchanged	-	Increased specificity of the 5'-primer, and exchanged control primer pair.
14	Added	Added	New primer pair for the B*5717 allele
16	Added	Added	New primer pair for the B*5718 allele

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*5701 to B*5719, i.e. all the currently recognized HLA-B*57 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-October-10, release 2.23.0.

RESOLUTION IN HLA-B*57 HOMO- AND HETEROZYGOTES:

The B*5701,5701 genotype give rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.