

Olerup SSP[®] HLA-B*57

Product number:	101.567-12 – including <i>Taq</i> polymerase 101.567-12u - without <i>Taq</i> polymerase
Lot number:	85M
Expiry date:	2014-May-01
Number of tests:	12
Number of wells per test:	24

CHANGES COMPARED TO THE PREVIOUS HLA-B*57 LOT (18L):

The HLA-B*57 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*57:01 to B*57:37, i.e. all the currently recognized HLA-B*57 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-B*57 HOMO- AND HETEROZYGOTES:

The B*57:01,57:01 genotype give rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.