

## **Olerup SSP<sup>®</sup> HLA-B\*57**

**Product number:** 101.567-12 – including *Taq* polymerase  
101.567-12u - without *Taq* polymerase  
**Lot number:** 62E  
**Expiry date:** 2010-May-01  
**Number of tests:** 12  
**Number of wells per test:** 16

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*57 LOT (31E):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
7	-	Modified	Increased specificity of specific primer pair.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

B\*5701 to B\*5716, i.e. all the currently recognized HLA-B\*57 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-April-08, release 2.21.0.

### **RESOLUTION IN HLA-B\*57 HOMO- AND HETEROZYGOTES:**

The B\*5701,5701 genotype give rise to a unique amplification pattern.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*57 SUBTYPINGS BY NON-HLA-B\*57 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.