

HLA-B*57
101.567-12 – licensed for PCR
101.567-12u – not licensed for PCR
Lot No.: Y36

1

Release Note - *Olerup* SSP™ HLA-B*57 – Lot No. Y36

Note -

Product number:	101.567-12 – licensed for PCR 101.567-12u - <u>not</u> licensed for PCR
Lot number:	Y36
Expiry date:	2009-September-01
Number of tests:	12
Number of tubes per test:	16

CHANGES COMPARED TO THE PREVIOUS HLA-B*57 LOT (X69):

Tube	5'-primer	3'-primer	rationale
11	-	Added	Primers added for the B*5707 and B*5714 alleles.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*5701 to B*5714, i.e. all the currently recognized HLA-B*57 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN HLA-B*57 HOMO- AND HETEROZYGOTES:

The B*5701,5701 genotype give rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.