

Olerup SSP® HLA-B*5701

Product number:	101.572-12 – including <i>Taq</i> polymerase 101.572-12u - without <i>Taq</i> polymerase
Lot number:	44G
Expiry date:	2011-October-01
Number of tests:	12
Number of wells per test:	8

CHANGES COMPARED TO THE PREVIOUS HLA-B*5701 Lot (88F):

The HLA-B*5701 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*570101 to B*570103, i.e. all the currently recognized HLA-B*5701 alleles, gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HLA-B*5701 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*5701 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.