

101.572-12 – including *Taq* polymerase101.572-12u – without *Taq* polymeraseLot No.: **88F****Olerup SSP[®] HLA-B*5701**

Product number: 101.572-12 – including *Taq* polymerase
101.572-12u - without *Taq* polymerase
Lot number: 88F
Expiry date: 2011-April-01
Number of tests: 12
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-B*5701 LOT (83E):

Well	5'-primer	3'-primer	rationale
3	-	Added	New 3'-primer for the B*5718 allele

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*570101 to B*570103, i.e. all the currently recognized HLA-B*5701 alleles, gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2009-January-16, release 2.24.0.

RESOLUTION IN HLA-B*5701 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*5701 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.